

ADHESION BOOSTERS AND INDIRECT BONDING

Enhance



ADHESION BOOSTER FOR ATYPICAL OR CONTAMINATED ENAMEL

Enhance will increase the adhesion of *Phase II*, *Rely a Bond*, *Light Bond* or *Excel* to fluorosed, hypo-calcified or deciduous enamel; acrylic composite or metal surface; any enamel that is close to the gingiva where minor surface moisture may reduce bond strength.

COMPONENTS

- 3cc Enhance Adhesion Booster part A (EBAA)
- 3cc Enhance Adhesion Booster part B (EBAB)
- 6cc Enhance Adhesion Booster part A (EBAA6)
- 6cc Enhance Adhesion Booster part B (EBAB6)

Therma Cure



HEAT CURABLE LAB ADHESIVE FOR CUSTOM BASE INDIRECT BONDING

Tacky Paste Viscosity. No bracket flotation. Unlimited Working Time. For bracket positioning and excess flash cleanup. *Heat Curable*. Place in toaster oven for 15 minutes at 325 degrees. *Highly Filled Composite*. Produces maximum bond strength with all metal, ceramic and plastic brackets. Simply bond in mouth with sealant.

Enhance L.C.



SINGLE COMPONENT LIGHT CURE

Provides all the benefits of chemical cure *Enhance*, but is optimized for use with *Light Bond Sealant*. This single liquid (no mixing) is applied to the etched enamel surface then lightly air dried. *Light Bond Sealant* is placed directly on the *Enhance L.C.* and light cured for 10 seconds (20 seconds is required if *Light Bond Filled Sealant* is used. Can only be used with *Light Bond Sealant*. Most convenient for *Light Bond* users.

Custom I.Q.



FLUORIDE RELEASING FAST SET INDIRECT BONDING SEALANT

Easy Application. Custom I.Q. Is an A & B sealant that is applied to the prepared tooth and custom bracket base without mixing. Simply apply part A to the custom bracket base and part B to the enamel and seat tray.

Fast Set Time. The indirect tray is held in position for 1 minute and then passively rests for 4 additional reducing the wait time in the custom indirect process.

Chemical Cure System. No curing light required. Bonds to any custom base pad when used with *Reliance Plastic Conditioner* or *Enhance Adhesion Booster*.

REFRIGERATION REQUIRED

Assure



UNIVERSAL BONDING RESIN

Assure Universal Bonding Resin will increase adhesion to normal, atypical, dry or slightly contaminated enamel. This universal bonding agent bonds to metal restorations without additional primers. No drying or light curing required and it can be used with any light or chemical cure system.