

Crown Lok

LUTING CEMENT FOR BANDS
 AND METAL CROWNS



Crown Lok is a self-curing two minute working time fluoride releasing cement. This versatile cement will bond chemically to all metal crowns, bands, and appliances, reducing the chance of failure and facilitating cleanup at the conclusion of treatment. Crown Lok is available in tooth shade or high visibility blue shade.

KIT CONTENTS

24gm paste A catalyst; 24gm paste B base; 9gm etching agent; mixing pad and spatulas.

COMPONENTS

- 6gm Crown Lok paste A in syringe (CLSA)
- 6gm Crown Lok paste B in syringe (CLSB)
- 6gm Crown Lok paste A in screw syringe (CLSAB)

These syringes mix and dispense paste in one easy step. The A & B syringes are joined and use a special mixing canula tip, which eliminates operator mixing and provides a consistent, automatic mix everytime.

KIT CONTENTS

Band-Lok (SMBLK) and Band-Lok Blue Shade (SMBLBK); 4 dual syringes of paste and 60 mixing tips.
 Crown Lok (SMCLBK): 4 dual paste syringes; 60 mixing tips and 9gm liquid etch.

Self Mixing Dual Cure Syringes

BAND-LOK, BAND-LOK BLUE, CROWN LOK



O.V.S.

ORTHODONTIC VENEER SYSTEM



O.V.S. is an excellent solution to a unique problem. O.V.S. has been perfected and used extensively by Dr. James Hilgers and Dr. Debra Cook in their own practice. A common problem when treating orthodontic cases is the need to bond an attachment to a peg lateral or an undersized tooth. With O.V.S. that problem is solved! O.V.S. will provide a complete in-office system that creates an esthetic and properly sized temporary veneer surface to bond your attachment. This saves you and your patients time and money by eliminating the need for a full temporary crown.

KIT CONTENTS

- Faci Form Crown Form #12 left (5 med and 5 large)
- Faci Form Crown Form #22 right (5 med and 5 large)
- (1) Ortho Veneer Paste Shade A-1
- (1) Ortho Veneer Paste Shade A-2
- (1) Ortho Veneer Paste Shade A-3
- (1) 3cc Light cure Sealant
- (20) Stick-n-Place instruments